# SID

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## References

1. Peters LB, Wesselink PR, Moorer WR. The fate and the role of bacteria left in root dentinal tubules. *Int Endod J* 1995; 28(2): 95-9.
2. Dalton BC, Orstavik D, Phillips C, Pettiette M, Trope M. Bacterial reduction with nickel-titanium rotary instrumentation. *J Endod* 1998; 24(11): 763-7.
3. Peters OA, Schonberger K, Laib A. Effects of four Ni-Ti preparation techniques on root canal geometry assessed by micro computed tomography. *Int Endod J* 2001; 34(3): 221-30.
4. Haapasalo M, Qian W. Irrigants and intracanal medicaments. In: Ingle JI, Bakland LK, editors. *Endodontics*. 6<sup>th</sup> ed. Hamilton, CA: BC Decker Inc; 2008. p. 996.
5. Hulsmann M, Hahn W. Complications during root canal irrigation--literature review and case reports. *Int Endod J* 2000; 33(3): 186-93.
6. Gomes BP, Souza SF, Ferraz CC, Teixeira FB, Zaia AA, Valdrighi L, et al. Effectiveness of 2% chlorhexidine gel and calcium hydroxide against *Enterococcus faecalis* in bovine root dentine in vitro. *Int Endod J* 2003; 36(4): 267-75.
7. Zerella JA, Fouad AF, Spangberg LS. Effectiveness of a calcium hydroxide and chlorhexidine digluconate mixture as disinfectant during retreatment of failed endodontic cases. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2005; 100(6): 756-61.
8. Manzur A, Gonzalez AM, Pozos A, Silva-Herzog D, Friedman S. Bacterial quantification in teeth with apical periodontitis related to instrumentation and different intracanal medications: a randomized clinical trial. *J Endod* 2007; 33(2): 114-8.
9. Saleh IM, Ruyter IE, Haapasalo M, Orstavik D. Survival of *Enterococcus faecalis* in infected dentinal tubules after root canal filling with different root canal sealers in vitro. *Int Endod J* 2004; 37(3): 193-8.
10. Barrette WC, Jr., Hannum DM, Wheeler WD, Hurst JK. General mechanism for the bacterial toxicity of hypochlorous acid: abolition of ATP production. *Biochemistry* 1989; 28(23): 9172-8.
11. Shahravan A, Haghdoust AA, Adl A, Rahimi H, Shadifar F. Effect of smear layer on sealing ability of canal obturation: a systematic review and meta-analysis. *J Endod* 2007; 33(2): 96-105.
12. McDonnell G, Russell AD. Antiseptics and disinfectants: activity, action, and resistance. *Clin Microbiol Rev* 1999; 12(1): 147-79.
13. McKenna SM, Davies KJ. The inhibition of bacterial growth by hypochlorous acid. Possible role in the bactericidal activity of phagocytes. *Biochem J* 1988; 254(3): 685-92.
14. Khademi A, Feizianfard M. The effect of EDTA and citric acid on smear layer removal of mesial canals of first mandibular molars, a scanning electron microscopic study. *J Res Med Sci* 2004; 9(2): 80-8.
15. Zehnder M, Kosicki D, Luder H, Sener B, Waltimo T. Tissue-dissolving capacity and antibacterial effect of buffered and unbuffered hypochlorite solutions. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2002; 94(6): 756-62.
16. Gomes BP, Ferraz CC, Vianna ME, Berber VB, Teixeira FB, Souza-Filho FJ. In vitro antimicrobial activity of several concentrations of sodium hypochlorite and chlorhexidine gluconate in the elimination of *Enterococcus faecalis*. *Int Endod J* 2001; 34(6): 424-8.
17. Radcliffe CE, Potouridou L, Qureshi R, Hababeh N, Qualtrough A, Worthington H, et al. Antimicrobial activity of varying concentrations of sodium hypochlorite on the endodontic microorganisms *Actinomyces israelii*, *A. naeslundii*, *Candida albicans* and *Enterococcus faecalis*. *Int Endod J* 2004; 37(7): 438-46.
18. Vianna ME, Gomes BP, Berber VB, Zaia AA, Ferraz CC, de Souza-Filho FJ. In vitro evaluation of the antimicrobial activity of chlorhexidine and sodium hypochlorite. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2004; 97(1): 79-84.
19. Waltimo TM, Orstavik D, Siren EK, Haapasalo MP. In vitro susceptibility of *Candida albicans* to four disinfectants and their combinations. *Int Endod J* 1999; 32(6): 421-9.
20. Haapasalo HK, Siren EK, Waltimo TM, Orstavik D, Haapasalo MP. Inactivation of local root canal medicaments by dentine: an in vitro study. *Int Endod J* 2000; 33(2): 126-31.
21. Portenier I, Waltimo T, Orstavik D, Haapasalo M. The susceptibility of starved, stationary phase, and growing cells of *Enterococcus faecalis* to endodontic medicaments. *J Endod* 2005; 31(5): 380-6.
22. Bystrom A, Sundqvist G. Bacteriologic evaluation of the effect of 0.5 percent sodium hypochlorite in endodontic therapy. *Oral Surg Oral Med Oral Pathol* 1983; 55(3): 307-12.
23. Bystrom A, Sundqvist G. The antibacterial action of sodium hypochlorite and EDTA in 60 cases of endodontic therapy. *Int Endod J* 1985; 18(1): 35-40.

24. Siqueira JF, Jr., Rocas IN, Santos SR, Lima KC, Magalhaes FA, de UM. Efficacy of instrumentation techniques and irrigation regimens in reducing the bacterial population within root canals. *J Endod* 2002; 28(3): 181-4.
25. Clegg MS, Vertucci FJ, Walker C, Belanger M, Britto LR. The effect of exposure to irrigant solutions on apical dentin biofilms in vitro. *J Endod* 2006; 32(5): 434-7.
26. Spangberg L, Engstrom B, Langeland K. Biologic effects of dental materials. 3. Toxicity and antimicrobial effect of endodontic antiseptics in vitro. *Oral Surg Oral Med Oral Pathol* 1973; 36(6): 856-71.
27. Sim TP, Knowles JC, Ng YL, Shelton J, Gulabivala K. Effect of sodium hypochlorite on mechanical properties of dentine and tooth surface strain. *Int Endod J* 2001; 34(2): 120-32.
28. Marending M, Luder HU, Brunner TJ, Knecht S, Stark WJ, Zehnder M. Effect of sodium hypochlorite on human root dentine--mechanical, chemical and structural evaluation. *Int Endod J* 2007; 40(10): 786-93.
29. Niu W, Yoshioka T, Kobayashi C, Suda H. A scanning electron microscopic study of dentinal erosion by final irrigation with EDTA and NaOCl solutions. *Int Endod J* 2002; 35(11): 934-9.
30. Baumgartner JC, Mader CL. A scanning electron microscopic evaluation of four root canal irrigation regimens. *J Endod* 1987; 13(4): 147-57.
31. Czonstkowsky M, Wilson EG, Holstein FA. The smear layer in endodontics. *Dent Clin North Am* 1990; 34(1): 13-25.
32. Baumgartner JC, Brown CM, Mader CL, Peters DD, Shulman JD. A scanning electron microscopic evaluation of root canal debridement using saline, sodium hypochlorite, and citric acid. *J Endod* 1984; 10(11): 525-31.
33. Loel DA. Use of acid cleanser in endodontic therapy. *J Am Dent Assoc* 1975; 90(1): 148-51.
34. Haapasalo M, Orstavik D. In vitro infection and disinfection of dentinal tubules. *J Dent Res* 1987; 66(8): 1375-9.
35. Orstavik D, Haapasalo M. Disinfection by endodontic irrigants and dressings of experimentally infected dentinal tubules. *Endod Dent Traumatol* 1990; 6(4): 142-9.
36. Haapasalo M, Shen Y, Qian W, Gao Y. Irrigation in endodontics. *Dent Clin North Am* 2010; 54(2): 291-312.
37. Machado-Silveiro LF, Gonzalez-Lopez S, Gonzalez-Rodriguez MP. Decalcification of root canal dentine by citric acid, EDTA and sodium citrate. *Int Endod J* 2004; 37(6): 365-9.
38. Gutmann JL, Saunders WP, Nguyen L, Guo IY, Saunders EM. Ultrasonic root-end preparation. Part 1. SEM analysis. *Int Endod J* 1994; 27(6): 318-24.
39. Hulsmann M, Heckendorff M, Lennon A. Chelating agents in root canal treatment: mode of action and indications for their use. *Int Endod J* 2003; 36(12): 810-30.
40. Zehnder M. Root canal irrigants. *J Endod* 2006; 32(5): 389-98.
41. Bjorvatn K, Skaug N, Selvig KA. Tetracycline-impregnated enamel and dentin: duration of antimicrobial capacity. *Scand J Dent Res* 1985; 93(3): 192-7.
42. Wikesjo UM, Baker PJ, Christersson LA, Genco RJ, Lyall RM, Hic S, et al. A biochemical approach to periodontal regeneration: tetracycline treatment conditions dentin surfaces. *J Periodontol* 1986; 21(4): 322-9.
43. Baker PJ, Evans RT, Coburn RA, Genco RJ. Tetracycline and its derivatives strongly bind to and are released from the tooth surface in active form. *J Periodontol* 1983; 54(10): 580-5.
44. Khademi AA, Mohammadi Z, Havaee A. Evaluation of the antibacterial substantivity of several intra-canal agents. *Aust Endod J* 2006; 32(3): 112-5.
45. Labahn R, Fahrenbach WH, Clark SM, Lie T, Adams DF. Root dentin morphology after different modes of citric acid and tetracycline hydrochloride conditioning. *J Periodontol* 1992; 63(4): 303-9.
46. Sae-Lim V, Wang CY, Choi GW, Trope M. The effect of systemic tetracycline on resorption of dried replanted dogs' teeth. *Endod Dent Traumatol* 1998; 14(3): 127-32.
47. Bridges JB, Owens PD, Stewart DJ. Tetracyclines and teeth. An experimental investigation into five types in the rat. *Br Dent J* 1969; 126(7): 306-11.
48. Russell AD. Activity of biocides against mycobacteria. *Journal of Applied Microbiology* 1996; 81(Supplement s25): 87S-101S.
49. Khademi A, Usefian E, Feizianfard M. Tissue dissolving ability of several endodontic irrigants on bovine pulp tissue. *Iran Endod J* 2007; 2(2): 65-8.
50. Russell AD, Day MJ. Antibacterial activity of chlorhexidine. *J Hosp Infect* 1993; 25(4): 229-38.
51. Heling I, Chandler NP. Antimicrobial effect of irrigant combinations within dentinal tubules. *Int Endod J* 1998; 31(1): 8-14.
52. Vahdaty A, Pitt Ford TR, Wilson RF. Efficacy of chlorhexidine in disinfecting dentinal tubules in vitro. *Endod Dent Traumatol* 1993; 9(6): 243-8.
53. Buck RA, Eleazer PD, Staat RH, Scheetz JP. Effectiveness of three endodontic irrigants at various tubular depths in human dentin. *J Endod* 2001; 27(3): 206-8.

54. Jeansonne MJ, White RR. A comparison of 2.0% chlorhexidine gluconate and 5.25% sodium hypochlorite as antimicrobial endodontic irrigants. *J Endod* 1994; 20(6): 276-8.
55. Zamany A, Safavi K, Spangberg LS. The effect of chlorhexidine as an endodontic disinfectant. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2003; 96(5): 578-81.
56. Ferraz CC, Gomes BP, Zaia AA, Teixeira FB, Souza-Filho FJ. In vitro assessment of the antimicrobial action and the mechanical ability of chlorhexidine gel as an endodontic irrigant. *J Endod* 2001; 27(7): 452-5.
57. Molander A, Reit C, Dahlen G. The antimicrobial effect of calcium hydroxide in root canals pretreated with 5% iodine potassium iodide. *Endod Dent Traumatol* 1999; 15(5): 205-9.
58. Ng YL, Spratt D, Sriskantharajah S, Gulabivala K. Evaluation of protocols for field decontamination before bacterial sampling of root canals for contemporary microbiology techniques. *J Endod* 2003; 29(5): 317-20.
59. Moller AJ. Microbiological examination of root canals and periapical tissues of human teeth. *Methodological studies. Odontol Tidskr* 1966; 74(5): Suppl-380.
60. Peciuliene V, Reynaud AH, Balciuniene I, Haapasalo M. Isolation of yeasts and enteric bacteria in root-filled teeth with chronic apical periodontitis. *Int Endod J* 2001; 34(6): 429-34.
61. Portenier I, Haapasalo H, Rye A, Waltimo T, Orstavik D, Haapasalo M. Inactivation of root canal medicaments by dentine, hydroxylapatite and bovine serum albumin. *Int Endod J* 2001; 34(3): 184-8.
62. Portenier I, Haapasalo H, Orstavik D, Yamauchi M, Haapasalo M. Inactivation of the antibacterial activity of iodine potassium iodide and chlorhexidine digluconate against *Enterococcus faecalis* by dentin, dentin matrix, type-I collagen, and heat-killed microbial whole cells. *J Endod* 2002; 28(9): 634-7.
63. Zehnder M, Schmidlin P, Sener B, Waltimo T. Chelation in root canal therapy reconsidered. *J Endod* 2005; 31(11): 817-20.
64. Marchesan MA, Pasternak JB, Afonso MM, Sousa-Neto MD, Paschoalato C. Chemical analysis of the flocculate formed by the association of sodium hypochlorite and chlorhexidine. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2007; 103(5): e103-e105.
65. Basrani BR, Manek S, Fillery E. Using diazotization to characterize the effect of heat or sodium hypochlorite on 2.0% chlorhexidine. *J Endod* 2009; 35(9): 1296-9.
66. Basrani BR, Manek S, Sodhi RN, Fillery E, Manzur A. Interaction between sodium hypochlorite and chlorhexidine gluconate. *J Endod* 2007; 33(8): 966-9.
67. Clarkson RM, Moule AJ, Podlich H, Kellaway R, Macfarlane R, Lewis D, et al. Dissolution of porcine incisor pulps in sodium hypochlorite solutions of varying compositions and concentrations. *Aust Dent J* 2006; 51(3): 245-51.
68. Williamson AE, Cardon JW, Drake DR. Antimicrobial susceptibility of monoculture biofilms of a clinical isolate of *Enterococcus faecalis*. *J Endod* 2009; 35(1): 95-7.
69. Shen Y, Qian W, Chung C, Olsen I, Haapasalo M. Evaluation of the effect of two chlorhexidine preparations on biofilm bacteria in vitro: a three-dimensional quantitative analysis. *J Endod* 2009; 35(7): 981-5.
70. Giardino L, Ambu E, Becce C, Rimondini L, Morra M. Surface tension comparison of four common root canal irrigants and two new irrigants containing antibiotic. *J Endod* 2006; 32(11): 1091-3.
71. Torabinejad M, Cho Y, Khademi AA, Bakland LK, Shabahang S. The effect of various concentrations of sodium hypochlorite on the ability of MTAD to remove the smear layer. *J Endod* 2003; 29(4): 233-9.
72. Torabinejad M, Khademi AA, Babagoli J, Cho Y, Johnson WB, Bozhilov K, et al. A new solution for the removal of the smear layer. *J Endod* 2003; 29(3): 170-5.
73. Shabahang S, Pouresmail M, Torabinejad M. In vitro antimicrobial efficacy of MTAD and sodium hypochlorite. *J Endod* 2003; 29(7): 450-2.
74. Shabahang S, Torabinejad M. Effect of MTAD on *Enterococcus faecalis*-contaminated root canals of extracted human teeth. *J Endod* 2003; 29(9): 576-9.
75. Johal S, Baumgartner JC, Marshall JG. Comparison of the antimicrobial efficacy of 1.3% NaOCl/BioPure MTAD to 5.25% NaOCl/15% EDTA for root canal irrigation. *J Endod* 2007; 33(1): 48-51.
76. Kho P, Baumgartner JC. A comparison of the antimicrobial efficacy of NaOCl/Biopure MTAD versus NaOCl/EDTA against *Enterococcus faecalis*. *J Endod* 2006; 32(7): 652-5.
77. Giardino L, Ambu E, Savoldi E, Rimondini R, Cassanelli C, Debbia EA. Comparative evaluation of antimicrobial efficacy of sodium hypochlorite, MTAD, and Tetraclean against *Enterococcus faecalis* biofilm. *J Endod* 2007; 33(7): 852-5.
78. Tay FR, Mazzoni A, Pashley DH, Day TE, Ngoh EC, Breschi L. Potential iatrogenic tetracycline staining of endodontically treated teeth via NaOCl/MTAD irrigation: a preliminary report. *J Endod* 2006; 32(4): 354-8.

## Irrigation solutions in endodontic treatment

Abbasali Khademi, Nastaran Farhadi, Ali Shokraneh\*

### Abstract

**Introduction:** *Elimination of pathogens from root canal system is the goal of endodontic treatment. Cleaning and shaping the root canal system is one of the most important steps in root canal therapy. Irrigation during root canal preparation leads to the elimination of microorganisms that cannot be eliminated by physical methods. In addition, the smear layer created during root canal preparation with hand and rotary instruments should be eliminated by irrigants. This review article summarizes the chemistry, biology and procedures for efficient and safe irrigation and presents information on recent developments in root canal irrigants. In addition, it evaluates irrigants from the viewpoint of their effect on microorganisms and the smear layer.*

**Materials and Methods:** *In this review article, studies were obtained and reviewed by running a search in various databases, including reference books, PubMed, ISI Web of Science, Scopus and Google search engine until 2011 using the terms root canal therapy, irrigant, antibacterial, chlorhexidine, smear layer and sodium hypochlorite.*

**Results:** *Several materials have been introduced as root canal irrigants. Although sodium hypochlorite is the most routine irrigant during root canal preparation, it has its limitations and disadvantages like other irrigants and is unable to eliminate the smear layer completely. Therefore, no irrigant can provide ideal properties of a suitable irrigant.*

**Conclusion:** *Use of two or more irrigants with an appropriate sequence has been recommended to eliminate microorganisms and the smear layer from the root canal system. Use of chlorhexidine as the final irrigant has substantivity properties.*

**Key words:** *Antibacterial agents, Chlorhexidine, Root canal therapy, Smear layer, Sodium hypochlorite*

**Received:** 7 Oct, 2012      **Accepted:** 8 Jan, 2013

**Address:** Postgraduate Student, Department of Endodontics, School of Dentistry, Isfahan University of Medical Sciences, Isfahan, Iran

**Email:** ali\_shokrane@dnt.mui.ac.ir

**Citation:** Khademi A, Farhadi N, Shokraneh A. **Irrigation solutions in endodontic treatment.** J Isfahan Dent Sch 2013; 8(6): 576-586.

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