



FLASH

ENDO POWER



📍 Bondent Group is a platform that provides integrated services for dental medical institutions.

We are committed both to improving world's level of dental medical technology and industry, and also to helping dentists succeed.

www.bondent.eu

✉ international@bondent.com

📞 **Europe:** +49 89 9545 5475

China: +86 187 5816 5369

📍 **Europe:** Bad-Berneck-Str.16
81549 München. Germany

China: 13th Floor, No. 655,
Fengzhou Road,
Jiading New Town,
Shanghai. China

📘 BondentEurope
BondentMEA
BondentIN
bondentTurkey


📷 @bondent_europe
@bondent_mea
@bondent_in
@bondent_turkey

Flash

Flash Story

Searching for a sufficiently safe reciprocation file, able to reduce harm to the natural tooth and during root canal preparation, Bondent Group started its studies on safe and excellent cutting efficiency reciprocation files. At the same time, Dr Nuno Pinto, from Portugal, joined us and gave support on design and clinical testing. We went through 2 years of a 5-time continuous sampling. in order to be submitted to laboratorial testing, which turned to result in ideal performance. We were assured to proceed with the making of clinical testing in cooperation with Dr Nuno Pinto. Due to the excellent clinical performance, Flash was able to be launched in 2022.

Clinical cases



Bondent KOL, Flash collaborator ,Nuno Pinto DDS,MSc


- Member of the Portuguese Society of Endodontology;
- Certified member of the European Society of Endodontology;
- Member of the American Association of Endodontists;
- Participant in the book "The root canal anatomy in permanent dentition";
- Professor of the Master's in Endodontics at Universitat Internacional de Catalunya.

"Endo power! It is a powerful file system with good flexibility and cutting efficiency. Due to the special flat design, it provides excellent debris removal." — Dr Nuno Pinto


Normal case




Easy case



Retreatment



Difficult case



Features

ISO color coding

Short

Short Shaft 11mm

Silicone stopper in the ISO color of the specific Flash instrument tip size



S cross section

Perfect cutting efficiency

Working length

Working length
21 mm, 25 mm and 31 mm

Flat side design

Flat side design increases the efficiency of debris removal


Heat treated NiTi material

Unique heat treatment technology: increases flexibility

Tip 6mm: no flat design, no blue nano coating

Minimally invasive

small size and taper design, in line with the concept of minimally invasive




Non-cutting tip



Product reference

Flash



20/04vt

21mm	UFLA004021020
25mm	UFLA004025020
31mm	UFLA004031020



25/06vt

21mm	UFLA006021025
25mm	UFLA006025025
31mm	UFLA006031025



35/05vt

21mm	UFLA005021035
25mm	UFLA005025035
31mm	UFLA005031035

Assortment


20/04vt, 25/06vt, 35/05vt

21mm	UFLA000021A01
25mm	UFLA000025A01
31mm	UFLA000031A01



iRoot Pro SRCT2000

Cordless Endo Motor
with integrated Apex Locator



Standard Gutta Percha Points

20/04	UDGP004000020
25/04	UDGP004000025
35/04	UDGP004000035