



# Small Device Big Clarity



**Bondream Air**

Portable Dental X-Ray Device

# INCREASED CONFIDENCE FOR PRECISE DIAGNOSIS

Focal Spot Size: 0.4 mm

🔍 Provide sharp and high-resolution images



# ADVANCED DC TECHNOLOGY FOR SAFER, SHARPER, SMARTER IMAGING



- 🔍 Delivering constant high voltage and current throughout every exposure, this advanced design provides stable and precise X-ray emissions



### One-Touch Operation for Maximum Efficiency

Pre-programmed exposure settings (e.g., for adult/child, anterior/posterior) allow fast, intuitive use with just a single button press—saving valuable clinical time

### User-Friendly Interface

Icon-based menu, intuitive layout make the interface easy to use—even for first-time users.

### Multiple Modes for Versatile Use

Available Voltage: 60kV/65kV/70kV

Available Current: 2.5mA/3.5mA

# INTUITIVE TOUCH SCREEN

# ONE TOUCH TO SELECT TOOTH POSITION

4 specialized imaging modes tailored to different tooth types, Anterior Teeth, Premolars, Molars, and Bitewing—ensuring optimized exposure settings for each region to deliver precise, high-quality diagnostic images.



## DUAL IMAGING COMPATIBILITY

Digital Sensor & Film Support in One Device



Technical Specification

| Model               | Bondream Air        |
|---------------------|---------------------|
| Weight              | 2.2 kg              |
| Focal spot          | 0.4mm               |
| Tube voltage        | 60 kV, 65 kV, 70 kV |
| Tube current        | 2.5mA, 3.5mA        |
| Exposure time range | 0.01-2.00S          |

Bondream Air





Distributed By.



**Bondent Group**  
Email: [international@bondent.com](mailto:international@bondent.com)



@Bondent Group